

## Uniblitz<sup>®</sup> FS25

25mm Uni-Stable Optical Shutter

### Overview

The Uniblitz FS series shutters are designed and optimized to operate directly from +5VDC and do not require a separate driver. Removing the +5VDC (0VDC) closes these shutters. A simple control circuit can be used to operate these shutter devices from a TTL trigger pulse. This control can be also accomplished with our new VLM1 TTL control interface (available soon).

This low-cost innovation provides the reliability of Uniblitz shutters (typical lifetime >300K cycles) at a single operating voltage.

### **Need Support?** Please <u>visit our website</u> or email us at <u>info@uniblitz.com</u>. Tel: <u>585-385-5930</u> | Toll-Free: <u>800-828-6972</u> | Fax: 585-385-6004 | 803 Linden Ave. Rochester, NY 14625

Updated 11/19 | Datasheet Version 5.1 | ©2019 Vincent Associates

### **Key Features**


- 25mm aperture
- Default closed operation, +5VDC opens the shutter, 0VDC closes the shutter
- RoHS Compliant
- Transfer time on opening:
   9.0 milliseconds
- Total opening time:
  - 16.0 milliseconds

## **Product Options**

FS25 2 3 4 5 6 - 7	Ex: FS25S2C0L-EC
-	<ul><li><b>2</b> Voltage:</li><li><b>5</b> Standard 5VDC</li></ul>
<ul> <li>1: Un-Housed</li> <li>2: Half-Housed</li> <li>3: Fully-Housed</li> </ul>	<ul> <li>Blades: <sup>1</sup></li> <li>C: Black carbon impregnated plastic (C-PET)</li> <li>T: Low Energy (Teflon<sup>®</sup>)</li> <li>ZM: High Energy (AlMgF2) <sup>2</sup></li> </ul>
• •	<ul><li>Connector:</li><li>L: 18" flying leads</li></ul>
<ul> <li>Encapsulated Coil:</li> <li>EC: Included</li> <li>Leave blank if not required</li> </ul>	

<sup>1</sup> Other blade coating options may be available by special order.
<sup>2</sup> Input side only; Teflon<sup>®</sup> coating is on opposite side to protect shutter blade surface. Light source must be input to the reflective side only.

# Shutter Timing



### FS25 (w/ 5VDC and C-PET blades)

### Time (msec.)


0 - A	Delay time on opening after current applied	7.0
A - C	Transfer time on opening	9.0
0 - C	Total opening time	16.0
C - E	Min. dwell time with min. input pulse	15.0
B - F	Min. equivalent exp. time	30.5
E - G	Transfer time on closing	22.0
A - G	Total window time	53.0
MET	Min. exposure time	30.0
TEP	Typical exposure pulse	>30.0

## **Technical Specifications**

Coil Resistance	Voltage to Open	Hold Voltage	
12 Ω	+5 VDC	+5 VDC	

<sup>1</sup> (Continuous/Burst) Continuous frequency rating specified at shutter's minimum exposure pulse. Burst frequency rating specified for four (4) seconds maximum with one (1) minute minimum between bursts.

Series	<b>Weight</b>	Operating	Max. Opening	Max. Closing	Max. Freq. of	Number of
	(Unhoused/Half/Housed)	Temp.	Bounce	Bounce	Operation <sup>1</sup>	Shutter Blades
FS25	22.0 g / 57.0 g / 92.5 g	-40 - +65 °C	15%	15%	5 Hz / 10 Hz	5




The synchronization system for FS shutter devices incorporates a small magnet mounted to the driving mechanism and a Hall effect sensor. When the device achieves approximately 80% of full open, the magnet causes the Hall effect sensor to change state, producing a signal to indicate that the shutter has switched to the active state. Shown to the left is the FS series shutter schematic which incorporates the electronic synchronization system. **There is no connection to the designated synchronization pins when an electronic sync. is not selected.** 


### Technical Drawings - FS25 (Half-Housed - FS25S2C0L)

#### NOTES:

- 1. UNITS: [MM] INCH.
- 2. LIMITED HIDDEN LINES SHOWN FOR CLARITY.
- 3. LEADS WIRES REMOVED FOR CLARITY.
- 4 Shutter actual aperture is 25.4mm, in mount aperture is reduced by female C-mount to 0.965 [24.5mm].



### Technical Drawings - FS25 (Un-Housed - FS25S1C0L)



## Technical Drawings - FS25 (Housed - FS25S3C0L)

#### NOTES:

- 1. UNITS: [MM] INCH
- 2. LIMITED HIDDEN LINES SHOWN FOR CLARITY.
- 3. SHUTTER APERTURE IS 25.4MM REDUCED TO (Ø.965 [24.5MM]) WHEN INSTALLED INTO HOUSING.
- 4. OPTIONAL REFLECTIVE SURFACE OPPOSITE ACTUATOR COIL SIDE.

